



WIPESPLUS TECHNICAL DATA

EFFICACY DATA FOR WIPESPLUS® HOSPITAL/SURGICAL WIPES EPA REG. NO. 70144-2-75399

PRODUCT DESCRIPTION

The unique formulation of this product is produced through a combination of various solvents, sequestering agents, chelating agents, nonionic detergents and other ingredients acting in synergy with three distinct antimicrobial active agents. This formulation may be safely applied directly to surfaces as an effective decontaminant cleaner and disinfectant intended for use in medical/surgical and animal care facilities.

Note: Do not use as a baby wipe.

CHEMICAL COMPOSITION Mycobactericidal /Bactericidal/Virucidal/Fungicidal/Gram Positive/Gram Negative

ACTIVE INGREDIENTS:

n-Alkyl (60% C14, 30% C16, 5% C12, 5% C18) dimethyl benzyl ammonium chloride.....	0.154%
n-Alkyl (68% C12, 32% C14) dimethyl ethylbenzyl ammonium chloride.....	0.154%
OTHER INGREDIENTS.....	99.692%
Total.....	100.00%

USE IN:

Patient Rooms/Emergency Rooms/Operating Rooms/Ambulatory Surgical Centers Neonatal & Infant Care Units/Child & Geriatric Care Facilities Plasma Pheresis & Blood Centers/Respiratory & Pulmonary Labs/Dialysis Facilities

FEATURES AND BENEFITS:

- The fast acting formula loosens and suspends blood and bodily fluids, non-organic soil and other cellular debris for easy cleaning.
- Contains no toxic chlorine or staining dyes and does not require special handling or storage.
- EPA registered WipesPlus® Hospital Surgical Wipes meet CDC & CMS guidelines for health care disinfectants and exceeds OSHA blood-borne pathogen requirements.
- Proven to kill Gram Positive and Gram Negative Bacteria, Tuberculosis, Chlamydia and Pathogenic Fungus.
- Effective against Surrogate Pathogens such as: Staphylococcus aureus, Salmonella (non-hospital) Pseudomonas (hospital), Trichopyton mentagrophytes, Herpes Simples, Rotavirus, Tuberculosis

CONTACT TIME:

Use a 2 minute contact time for disinfection against all virus, bacteria, and fungi claimed.

EFFECTIVE AGAINST:

- Acinetobacter baumannii (ATCC 19606)
- Clostridium difficile (Veg) (ATCC 9689)
- Escherichia coli (E. Coli ESBL) (ATCC BAA-196)
- Hepatitis B Virus (HBV)
- Hepatitis C Virus (HCV)
- Herpes Simplex Virus (HSV)
- Human Immunodeficiency Virus (HIV-1)
- Influenza A.Virus
- Klebsiella pneumoniae –Carbapenem Resistant Strain (KPC)
- Listeria (ATCC 7644)
- Methicillin Resistant Staphylococcus Aureus (MRSA) (ATCC 33591)
- Mycobacterium tuberculosis (bovis) (TB)
- Pseudomonas aeruginosa
- Rhinovirus
- Rotavirus
- Salmonella entérica
- Serratia marcescens
- Staphylococcus aureus
- Swine Influenza A Virus (H1N1 strain)
- Trichophyton mentagrophytes
- Vancomycin Resistant Enterococcus (VRE)

QUANTITY:

Item No. 33906 240ct 12 per case 6"x 6.75"

SAFETY (See product for full safety precautions):

Caution: Harmful if swallowed. Harmful if inhaled.

Avoid breathing vapors. Causes moderate eye irritation.