


# GLOVE SIZING GUIDE

## Method 1

Using a measuring tape, measure the circumference of your hand just below the knuckles. Use the chart below to determine your size

| Measurement | Glove Size |   |
|-------------|------------|---|
| 5"          | XXXS       | ● |
| 6"          | XXS        | ● |
| 7"          | XS         | ● |
| 8"          | S          | ● |
| 9"          | M          | ● |
| 10"         | L          | ● |
| 11"         | XL         | ● |
| 12"         | XXL        | ● |
| 13"         | XXXL       | ● |



## Method 2

Print this page making sure that it is at 100% (NOT scaled to fit).

Place your right hand on the scale below with the side of your index finger to the right of the red line. The number closest to the right edge of your hand will be the recommended glove size.

