

**974365 Compatibility Notes - 02/26/2016**

<b>Model</b>	
T-49F	Standard after 1/13/14 per pcn 7415. Freezer does not require standoff. Will fit on older cabinets, but it is not recommended to mix with white shelves.
TS-49F	Standard after 1/13/14 per pcn 7415. Freezer does not require standoff. Will fit on older cabinets, but it is not recommended to mix with white shelves.
T-49	Standard approximately after 11/2/14 per pcn 7955. Will fit on older cabinets, but it is not recommended to mix with white shelves.
TS-49G	Standard after 1/13/14 per pcn 7415. Will fit on older cabinets, but it is not recommended to mix with white shelves.
T-49G	Standard approximately after 11/2/14 per pcn 7955. Will fit on older cabinets, but it is not recommended to mix with white shelves.
T-49-2-G-2	Standard approximately after 11/2/14 per pcn 7955. Will fit on older cabinets, but it is not recommended to mix with white shelves.
TS-49	Standard after 1/13/14 per pcn 7415. Will fit on older cabinets, but it is not recommended to mix with white shelves.
T-49-4	Standard approximately after 11/2/14 per pcn 7955. Will fit on older cabinets, but it is not recommended to mix with white shelves.
T-49-HC	Standard approximately after 11/2/14 per pcn 7955. Will fit on older cabinets, but it is not recommended to mix with white shelves.
T-49G-HC-LD	Standard approximately after 11/2/14 per pcn 7955.
T-49F-HC	Standard approximately after 11/2/14 per pcn 7955. Will fit on older cabinets, but it is not recommended to mix with white shelves.