



Conversion of Active Dry Yeast to SAF Red (15909, 31105, 31150)

Active Dry Yeast (ounces)	SAF Red Instant Dry Yeast (ounces)	Water Adjustment (ounces)
8	8	0
16	16	0
Conversion	1 to 1	

Conversion Example: When converting 8 ounces of Red Star active dry to SAF Red instant dry yeast. SAF Red instant dry yeast = 8 ounces x 1 = 8 ounces. You do not need to adjust the amount of water you are converting.

Active Dry Yeast (grams)	SAF Red Instant Dry Yeast (grams)	Water Adjustment (grams)
100	100	0
250	250	0
Conversion	1 to 1	

Conversion Example: When converting 100 grams of Red Star active dry to SAF Red instant dry yeast. SAF Red instant dry yeast = 100 grams x 1 = 100 grams. You do not need to adjust the amount of water you are converting.

- * Usage level may vary depending on application. The above is a suggested guideline only *
- * Store any unused SAF Red instant dry yeast in airtight container preferably in refrigeration *
- * If rehydrating in water, suggested water temperature is between 90 – 100 °F. *

Conversión de Levadura Seca Activa a SAF Roja (15909, 31105, 31150)

Levadura Seca Activa (onzas)	SAF Roja Levadura Seca Instantánea (onzas)	Ajuste de Agua (onzas)
8	8	0
16	16	0
Conversión	1 a 1	

Ejemplo de Conversión: Para convertir 8 onzas de levadura activa Red Star a SAF Roja levadura instantánea. SAF Roja levadura instantánea = 8 onzas x 1 = 8 onzas. No es necesario ajustar la cantidad de agua que usted está convirtiendo.

Levadura Seca Activa (gramos)	SAF Roja Levadura Seca Instantánea (gramos)	Ajuste de Agua (gramos)
100	100	0
250	250	0
Conversión	1 a 1	

Ejemplo de Conversión: Para convertir 100 gramos de levadura activa Red Star a SAF Roja levadura instantánea. SAF Roja levadura instantánea = 100 gramos x 1 = 100 gramos. No es necesario ajustar la cantidad de agua que usted está convirtiendo.

- * El nivel de uso puede variar según la aplicación. Lo anterior es solo una guía sugerida *
- * Guarde cualquier SAF Roja instantánea sin utilizar en recipiente a prueba de aire en refrigeración *
- * Si va a rehidratar en agua, temperatura sugerida del agua es de 32 – 38 °C. *

For technical support, please call: 414-615-4094 or 414-213-2732 (9:00 am – 5:00 pm CDT)